**KERAMIT 450 +**

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**INSTRUCTIONS FOR USE**

**Spruing** - Sprues with a diameter of at least 3 mm are recommended. As feeding technique both direct and indirect techniques can be used: in case of bar casting technique, the diameter of the bar should be 4 mm and 2 mm for the feeder sprues. Be sure that the model is out of the heat center. The minimum wall thickness of the model should be of 0,4 mm.

**Investment** - Use a phosphate bonded investment for ceramic alloys, carefully following the manufacturer's instructions. Final temperature of pre-heating should be of 820°C, to be held during 30-60 minutes, depending on the dimensions and the number of the cylinders.

**Melting** - Various melting techniques can be used, always adding at least 50% of new alloy. It is suggested to use ceramic crucibles, without using flux. When melting with a torch, use the neutral zone (reducing) of the oxygen-propane flame; the alloy is ready to cast when the ingots pool and roll. Follow the same procedure when melting with an electric or induction centrifuge.

**Casting temperature:** 1400°C.

Allow to bench cool to room temperature.

**Finishing and cleaning** - Use abrasive instruments made of ceramic binder. Sand blast the metal with new aluminium oxide (50 microns). Clean with steam or with an ultrasonic machine for 10 minutes.

**Metal conditioning** - Oxidation is to be made in furnace at 650°C increasing the temperature till 980°C, in air, and keeping for 10 minutes.

**Opaque** - Apply opaque directly on the oxidized layer of the framework. A first layer very liquid (wash brand) and fired at a temperature kept 10°C over the usual one is suggested. Then proceed as usual.

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